

## DNJB2 Polyclonal Antibody

Cat No: HR1AP11793

For research use only

### Overview

Product Name	DNJB2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	DNAJB2 HSJ1 HSPF3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DnaJ heat shock protein family (Hsp40) member B2(DNAJB2) Homo sapiens This gene is almost exclusively expressed in the brain, mainly in the neuronal layers. It encodes a protein that shows sequence similarity to bacterial DnaJ protein and the yeast homologs. In bacteria, this protein is implicated in protein folding and protein complex dissociation. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2011],