

## CCBE1 Polyclonal Antibody

Cat No: HR1AP11787

For research use only

### Overview

Product Name	CCBE1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	CCBE1 KIAA1983 UNQ1921/PRO4395
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	collagen and calcium binding EGF domains 1(CCBE1) Homo sapiens This gene is thought to function in extracellular matrix remodeling and migration. It is predominantly expressed in the ovary, but down regulated in ovarian cancer cell lines and primary carcinomas, suggesting its role as a tumour suppressor. Mutations in this gene have been associated with Hennekam lymphangiectasia-lymphedema syndrome, a generalized lymphatic dysplasia in humans. [provided by RefSeq, Mar 2010],