

T184B Polyclonal Antibody

Cat No: HR1AP11782

For research use only

Overview

Product Name	T184B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	TMEM184B C22orf5 PSEC0108
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May activate the MAP kinase signaling pathway.,similarity:Belongs to the UPF0206 family.,