

## RTDR1 Polyclonal Antibody

Cat No: HR1AP11778

For research use only

### Overview

Product Name	RTDR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	RTDR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	radial spoke head 14 homolog(RSPH14) Homo sapiens This gene encodes a protein with no known function but with slight similarity to a yeast vacuolar protein. The gene is located in a region deleted in pediatric rhabdoid tumors of the brain, kidney and soft tissues, but mutations in this gene have not been associated with the disease. [provided by RefSeq, Jul 2008],