

## TNF18 Polyclonal Antibody

Cat No: HR1AP11780

For research use only

### Overview

Product Name	TNF18 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	TNFSF18 AITRL GITRL TL6 UNQ149/PRO175
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tumor necrosis factor superfamily member 18(TNFSF18) Homo sapiens The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptor TNFRSF18/AITR/GITR. It has been shown to modulate T lymphocyte survival in peripheral tissues. This cytokine is also found to be expressed in endothelial cells, and is thought to be important for interaction between T lymphocytes and endothelial cells. [provided by RefSeq, Jul 2008],