

## STAP2 Polyclonal Antibody

Cat No: HR1AP11776

For research use only

### Overview

Product Name	STAP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	STAP2 BKS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>signal transducing adaptor family member 2(STAP2) Homo sapiens This gene encodes the substrate of breast tumor kinase, an Src-type non-receptor tyrosine kinase. The encoded protein possesses domains and several tyrosine phosphorylation sites characteristic of adaptor proteins that mediate the interactions linking proteins involved in signal transduction pathways. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008],</p>