

BIN2 Polyclonal Antibody

Cat No: HR1AP11775

For research use only

Overview

Product Name	BIN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	62kDa
GeneID?Human?	BIN2 BRAP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	induction:Up-regulated during monocytic differentiation.,similarity:Contains 1 BAR domain.,subunit:Interacts with BIN1.,tissue specificity:Highest level expression seen in spleen and peripheral blood leukocytes and is also expressed at high levels in thymus, colon and placenta, suggesting preferential expression in hematopoietic tissues.,