

## TF2L1 Polyclonal Antibody

Cat No: HR1AP11771

For research use only

### Overview

Product Name	TF2L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	TFCP2L1 CRTR1 LBP9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Transcriptional suppressor. May suppress UBP1-mediated transcriptional activation. Modulates the placental expression of CYP11A1.,similarity:Belongs to the grh/CP2 family. CP2 subfamily.,tissue specificity:Highly expressed in placental JEG-3 cells and very low levels of expression in non-steroidogenic cells. No expression was seen in adrenal NCI-H295A cells or in adrenal tissue.,