

TOPRS Polyclonal Antibody

Cat No: HR1AP11769

For research use only

Overview

Product Name	TOPRS Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	114kDa
GeneID?Human?	TOPORS LUN TP53BPL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>TOP1 binding arginine/serine rich protein(TOPORS) Homo sapiens This gene encodes a nuclear protein which is serine and arginine rich, and contains a RING-type zinc finger domain. It is highly expressed in the testis, and functions as an ubiquitin-protein E3 ligase. Mutations in this gene are associated with retinitis pigmentosa type 31. Alternatively spliced transcript variants, encoding different isoforms, have been observed for this locus. [provided by RefSeq, Sep 2010],</p>