

ZNF71 Polyclonal Antibody

Cat No: HR1AP11768

For research use only

Overview

Product Name	ZNF71 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	ZNF71 EZFIT
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in transcriptional regulation.,induction:By tumor necrosis factor alpha (TNFalpha).,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 13 C2H2-type zinc fingers.,tissue specificity:Highly expressed in placenta, followed by brain, testis, pancreas, heart, small intestine, muscle, uterus, prostate and peripheral blood leukocytes. Not detected in liver, lung, colon, stomach, salivary and thyroid gland.,