

ESRP2 Polyclonal Antibody

Cat No: HR1AP11763

For research use only

Overview

Product Name	ESRP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	ESRP2 RBM35B PP7059
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	epithelial splicing regulatory protein 2(ESRP2) Homo sapiens ESRP2 is an epithelial cell-type-specific splicing regulator (Warzecha et al., 2009 [PubMed 19285943]).[supplied by OMIM, Aug 2009],