

# MOB1A Polyclonal Antibody

Cat No: HR1AP11765

For research use only

## Overview

Product Name	MOB1A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	MOB1A C2orf6 MOB4B MOBK1B MOBKL1B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>MOB kinase activator 1A(MOB1A) Homo sapiens The protein encoded by this gene is a component of the Hippo signaling pathway, which controls organ size and tumor growth by enhancing apoptosis. Loss of the encoded protein results in cell proliferation and cancer formation. The encoded protein is also involved in the control of microtubule stability during cytokinesis. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2015],</p>