

MANS1 Polyclonal Antibody

Cat No: HR1AP11764

For research use only

Overview

Product Name	MANS1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	MANSC1 LOH12CR3 UNQ316/PRO361
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 MANSC domain.,tissue specificity:Widely expressed.,