

TSYL2 Polyclonal Antibody

Cat No: HR1AP11759

For research use only

Overview

Product Name	TSYL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	TSPYL2 CDA1 DENTT TSPX HRIHFB2216
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TSPY like 2(TSPYL2) Homo sapiens This gene encodes a member of the testis-specific protein Y-encoded, TSPY-like/SET/nucleosome assembly protein-1 superfamily. The encoded protein is localized to the nucleolus where it functions in chromatin remodeling and as an inhibitor of cell-cycle progression. This protein may play a role in the suppression of tumor growth. [provided by RefSeq, Sep 2009],