

BACH2 Polyclonal Antibody

Cat No: HR1AP11754

For research use only

Overview

Product Name	BACH2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	BACH2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Transcriptional regulator that acts as repressor or activator. Binds to Maf recognition elements (MARE). Play important roles in coordinating transcription activation and repression by MAFK.,similarity:Belongs to the bZIP family. CNC subfamily.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 1 bZIP domain.,subunit:Heterodimer of BACH2 and Maf-related transcription factors.,tissue specificity:B-cell specific.,