

TRI31 Polyclonal Antibody

Cat No: HR1AP11756

For research use only

Overview

Product Name	TRI31 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	TRIM31 C6orf13 RNF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tripartite motif containing 31(TRIM31) Homo sapiens This gene encodes a protein that functions as an E3 ubiquitin-protein ligase. This gene shows altered expression in certain tumors and may be a negative regulator of cell growth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],