

T106C Polyclonal Antibody

Cat No: HR1AP11752

For research use only

Overview

Product Name	T106C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	TMEM106C EMOG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the TMEM106 family.,