

TCHP Polyclonal Antibody

Cat No: HR1AP11750

For research use only

Overview

Product Name	TCHP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	TCHP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May act as a 'capping' or 'branching' protein for keratin filaments in the cell periphery. May regulate K8/K18 filament and desmosome organization mainly at the apical or peripheral regions of simple epithelial cells.,similarity:Belongs to the TCHP family.,subunit:Interacts specifically with keratin proteins including, KRT5, KRT6A, KRT8, KRT14, KRT16 and KRT18.,