

# EBP2 Polyclonal Antibody

Cat No: HR1AP11747

For research use only

## Overview

|                       |  |
|-----------------------|--|
| Product Name          | EBP2 Polyclonal Antibody   |
| Source                | Rabbit   |
| Applications          | WB,ELISA   |
| Species Reactivity    | Human,Mouse  |
| Recommended Dilutions |  |
| Immunogen             |  |
| Species               | Rabbit   |
| Storage               | -20°C/1 year   |
| Isotype               |  |
| Clonality             |  |
| Concentration         | 1 mg/ml  |
| Observed band         | 33kDa  |
| GeneID?Human?         | EBNA1BP2 EBP2  |
| Human Swiss-Prot No.  |  |
| Cellular localization |  |
| Alternative Names     |  |
| Background            | function:Required for the processing of the 27S pre-rRNA.,similarity:Belongs to the EBP2 family.,subunit:Specifically interacts with EBV EBNA1. The EBNA1-EBP2 interaction is important for the stable segregation of EBV episomes during cell division.,tissue specificity:Ubiquitous., |