

DIRC2 Polyclonal Antibody

Cat No: HR1AP11744

For research use only

Overview

Product Name	DIRC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	DIRC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>disrupted in renal carcinoma 2(DIRC2) Homo sapiens This gene encodes a membrane-bound protein from the major facilitator superfamily of transporters. Disruption of this gene by translocation has been associated with haplo-insufficiency and renal cell carcinomas. Alternatively spliced transcript variants have been described, but their biological validity has not yet been determined. [provided by RefSeq, Jul 2008],</p>