

NEB2 Polyclonal Antibody

Cat No: HR1AP11743

For research use only

Overview

Product Name	NEB2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	89kDa
GeneID?Human?	PPP1R9B PPP1R6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>protein phosphatase 1 regulatory subunit 9B(PPP1R9B) Homo sapiens This gene encodes a scaffold protein that functions as a regulatory subunit of protein phosphatase 1a. Expression of this gene is particularly high in dendritic spines, suggesting that the encoded protein may play a role in receiving signals from the central nervous system. The encoded protein has putative tumor suppressor function and decreased expression has been observed in tumors. [provided by RefSeq, Feb 2014],</p>