

ALKB3 Polyclonal Antibody

Cat No: HR1AP11740

For research use only

Overview

Product Name	ALKB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	ALKBH3 ABH3 DEPC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>alkB homolog 3, alpha-ketoglutaratedependent dioxygenase(ALKBH3) Homo sapiens The Escherichia coli AlkB protein protects against the cytotoxicity of methylating agents by repair of the specific DNA lesions generated in single-stranded DNA. ALKBH2 (MIM 610602) and ALKBH3 are E. coli AlkB homologs that catalyze the removal of 1-methyladenine and 3-methylcytosine (Duncan et al., 2002 [PubMed 12486230]).[supplied by OMIM, Mar 2008],</p>