



## DUX5 Polyclonal Antibody

Cat No: HR1AP11739

For research use only

### Overview

Product Name	DUX5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	DUX5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	double homeobox 5(DUX5) Homo sapiens The human genome contains hundreds of repeats of the 3.3-kb family in regions associated with heterochromatin. The DUX gene family, including DUX5, resides within these 3.3-kb repeated elements (Beckers et al., 2001 [PubMed 11245978]). See DUX4 (MIM 606009).[supplied by OMIM, Mar 2008],