

## DMRTB Polyclonal Antibody

Cat No: HR1AP11736

For research use only

### Overview

Product Name	DMRTB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	DMRTB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue specificity:Testis.,