

KAD8 Polyclonal Antibody

Cat No: HR1AP11737

For research use only

Overview

Product Name	KAD8 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	AK8 C9orf98
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:ATP + AMP = 2 ADP.,similarity:Belongs to the adenylate kinase family.,