

## E2F7 Polyclonal Antibody

Cat No: HR1AP11734

For research use only

### Overview

Product Name	E2F7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	E2F7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	E2F transcription factor 7(E2F7) Homo sapiens E2F transcription factors, such as E2F7, play an essential role in the regulation of cell cycle progression (Di Stefano et al., 2003 [PubMed 14633988]).[supplied by OMIM, May 2008],