



# MINA Polyclonal Antibody

Cat No: HR1AP11735

For research use only

## Overview

Product Name	MINA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	MINA MDIG MINA53 NO52
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MYC induced nuclear antigen(MINA) Homo sapiens MINA is a c-Myc (MYC; MIM 190080) target gene that may play a role in cell proliferation or regulation of cell growth. (Tsuneoka et al., 2002 [PubMed 12091391]; Zhang et al., 2005 [PubMed 15897898]).[supplied by OMIM, May 2008],