

ARL11 Polyclonal Antibody

Cat No: HR1AP11733

For research use only

Overview

Product Name	ARL11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	ARL11 ARLTS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ADP ribosylation factor like GTPase 11(ARL11) Homo sapiens This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers. [provided by RefSeq, May 2010],