

DIRC1 Polyclonal Antibody

Cat No: HR1AP11732

For research use only

Overview

Product Name	DIRC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	DIRC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	disrupted in renal carcinoma 1(DIRC1) Homo sapiens This gene consists of two exons and encodes a 104aa protein. Disruption of this gene by translocation t(2;3)(q33;q21) is associated with familial clear cell renal cancer. [provided by RefSeq, Sep 2010],