

## SNPC3 Polyclonal Antibody

Cat No: HR1AP11730

For research use only

### Overview

Product Name	SNPC3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	SNAPC3 SNAP50
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,subunit:Part of the SNAPc complex composed of 5 subunits: SNAPC1, SNAPC2, SNAPC3, SNAPC4 and SNAPC5. SNAPC3 interacts with SNAPC1.,</p>