

ZN713 Polyclonal Antibody

Cat No: HR1AP11726

For research use only

Overview

Product Name	ZN713 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	ZNF713
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger protein 713(ZNF713) Homo sapiens The protein encoded by this gene contains C2H2 zinc finger domains. In some individuals, a CGG-repeat expansion from 5-22 repeats to 68-450 repeats has been identified in the first intron of this gene. This mutation is thought to effect the expression of this gene and it has been proposed that it may be associated with Autistic Spectrum Disorder. [provided by RefSeq, Jul 2016],