

CPXCR Polyclonal Antibody

Cat No: HR1AP11723

For research use only

Overview

Product Name	CPXCR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	CPXCR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>CPX chromosome region, candidate 1(CPXCR1) Homo sapiens This gene is one of several genes identified in a region of the X chromosome associated with an X-linked cleft palate (CPX) disorder. The encoded protein contains a motif similar to a motif found in zinc-finger proteins. Mutation analysis of this gene has not revealed any mutation which causes the CPX disorder. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Sep 2011].</p>