

# MKX Polyclonal Antibody

Cat No: HR1AP11717

For research use only

## Overview

Product Name	MKX Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	MKX C10orf48 IRXL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mohawk homeobox(MKX) Homo sapiens The protein encoded by this gene is an IRX family-related homeobox protein that may play a role in cell adhesion. Studies in mice suggest that this protein may be a regulator of tendon development. Two transcript variants encoding the same protein have been found for this gene.[provided by RefSeq, Jun 2011],