



HOMEZ Polyclonal Antibody

Cat No: HR1AP11714

For research use only

Overview

Product Name	HOMEZ Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	HOMEZ KIAA1443
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May function as a transcriptional regulator.,similarity:Contains 2 homeobox DNA-binding domains.,subunit:Homodimer or heterodimer (Potential). Interacts with HOXC8.,tissue specificity:Ubiquitous. Strongly expressed in adult testis and kidney as well as fetal lung and kidney.,