

AEBP1 Polyclonal Antibody

Cat No: HR1AP11712

For research use only

Overview

Product Name	AEBP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	127kDa
GeneID?Human?	AEBP1 ACLP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>AE binding protein 1(AEBP1) Homo sapiens This gene encodes a member of carboxypeptidase A protein family. The encoded protein may function as a transcriptional repressor and play a role in adipogenesis and smooth muscle cell differentiation. Studies in mice suggest that this gene functions in wound healing and abdominal wall development. Overexpression of this gene is associated with glioblastoma. [provided by RefSeq, May 2013],</p>