

GKN2 Polyclonal Antibody

Cat No: HR1AP11710

For research use only

Overview

Product Name	GKN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	GKN2 BLOT GDDR TFIZ1 UNQ465/PRO813
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>gastrokine 2(GKN2) Homo sapiens The secretory protein encoded by this gene is produced in gastric surface mucous cells, where it can bind trefoil factor family peptide 1 or gastrokine-1. This gene may be a tumor suppressor gene, as its expression is markedly decreased in gastric cancer tissues. The encoded protein interacts with gastrokine-1 and regulates homeostasis of the gastric mucosa. [provided by RefSeq, Dec 2015].</p>