

RB6I2 Polyclonal Antibody

Cat No: HR1AP11711

For research use only

Overview

Product Name	RB6I2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	122kDa
GeneID?Human?	ERC1 ELKS KIAA1081 RAB6IP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ELKS/RAB6-interacting/CAST family member 1(ERC1) Homo sapiens The protein encoded by this gene is a member of a family of RIM-binding proteins. RIMs are active zone proteins that regulate neurotransmitter release. This gene has been found fused to the receptor-type tyrosine kinase gene RET by gene rearrangement due to the translocation t(10;12)(q11;p13) in thyroid papillary carcinoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014],</p>