

## POGZ Polyclonal Antibody

Cat No: HR1AP11705

For research use only

### Overview

Product Name	POGZ Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	155kDa
GeneID?Human?	POGZ KIAA0461 SUHW5 ZNF280E ZNF635 Nbla00003
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	pogo transposable element with ZNF domain(POGZ) Homo sapiens The protein encoded by this gene appears to be a zinc finger protein containing a transposase domain at the C-terminus. This protein was found to interact with the transcription factor SP1 in a yeast two-hybrid system. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Aug 2010],