

# LAR1B Polyclonal Antibody

Cat No: HR1AP11697

For research use only

## Overview

Product Name	LAR1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	LARP1B LARP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>La ribonucleoprotein domain family member 1B(LARP1B) Homo sapiens This gene encodes a protein containing domains found in the La related protein of Drosophila melanogaster. La motif-containing proteins are thought to be RNA-binding proteins, where the La motif and adjacent amino acids fold into an RNA recognition motif. The La motif is also found in proteins unrelated to the La protein. Alternative splicing has been observed at this locus and multiple variants, encoding distinct isoforms, are described. Additional splice variation has been identified but the full-length nature of these transcripts has not been determined. [provided by RefSeq, Jun 2013].</p>