

DMRTA Polyclonal Antibody

Cat No: HR1AP11694

For research use only

Overview

Product Name	DMRTA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	DMRTA1 DMO
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in sexual development.,similarity:Belongs to the DMRT family.,similarity:Contains 1 DM DNA-binding domain.,tissue specificity:Expressed in liver, kidney, pancreas, prostate and weakly detected in testis and ovary.,