

## HES5 Polyclonal Antibody

Cat No: HR1AP11692

For research use only

### Overview

Product Name	HES5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	HES5 BHLHB38
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	hes family bHLH transcription factor 5(HES5) Homo sapiens This gene encodes a member of a family of basic helix-loop-helix transcriptional repressors. The protein product of this gene, which is activated downstream of the Notch pathway, regulates cell differentiation in multiple tissues. Disruptions in the normal expression of this gene have been associated with developmental diseases and cancer. [provided by RefSeq, Dec 2008],