

SNPC4 Polyclonal Antibody

Cat No: HR1AP11691

For research use only

Overview

Product Name	SNPC4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	161kDa
GeneID?Human?	SNAPC4 SNAP190
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>small nuclear RNA activating complex polypeptide 4(SNAPC4) Homo sapiens This gene encodes the largest subunit of the small nuclear RNA-activating protein (SNAP) complex. The encoded protein contains a Myb DNA-binding domain, and is essential for RNA polymerase II and III polymerase transcription from small nuclear RNA promoters. A mutation in this gene is associated with ankylosing spondylitis. [provided by RefSeq, Jul 2016].</p>