

GSAML Polyclonal Antibody

Cat No: HR1AP11690

For research use only

Overview

Product Name	GSAML Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	GCSAML C1orf150
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	germinal center associated signaling and motility like(GCSAML) Homo sapiens This gene encodes a protein thought to be a signaling molecule associated with germinal centers, the sites of proliferation and differentiation of mature B lymphocytes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013],