

TAF7L Polyclonal Antibody

Cat No: HR1AP11688

For research use only

Overview

Product Name	TAF7L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	TAF7L TAF2Q
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TATA-box binding protein associated factor 7 like(TAF7L) Homo sapiens This gene is similar to a mouse gene that encodes a TATA box binding protein-associated factor, and shows testis-specific expression. The encoded protein could be a spermatogenesis-specific component of the DNA-binding general transcription factor complex TFIID. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009],