

TAF1B Polyclonal Antibody

Cat No: HR1AP11687

For research use only

Overview

Product Name	TAF1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	TAF1B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TATA-box binding protein associated factor, RNA polymerase I subunit B(TAF1B) Homo sapiens Initiation of transcription by RNA polymerase I requires the formation of a complex composed of the TATA-binding protein (TBP) and three TBP-associated factors (TAFs) specific for RNA polymerase I. This complex, known as SL1, binds to the core promoter of ribosomal RNA genes to position the polymerase properly and acts as a channel for regulatory signals. This gene encodes one of the SL1-specific TAFs. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016].