



ZN750 Polyclonal Antibody

Cat No: HR1AP11684

For research use only

Overview

Product Name	ZN750 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	ZNF750
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger protein 750(ZNF750) Homo sapiens This gene encodes a protein with a nuclear localization site and a C2H2 zinc finger domain. Mutations in this gene have been associated with seborrhea-like dermatitis with psoriasiform elements. [provided by RefSeq, Jul 2008],