

## HEN2 Polyclonal Antibody

Cat No: HR1AP11682

For research use only

### Overview

Product Name	HEN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	NHLH2 BHLHA34 HEN2 KIAA0490
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May serve as DNA-binding protein and may be involved in the control of cell-type determination, possibly within the developing nervous system.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,