

ADCY3 Polyclonal Antibody

Cat No: HR1AP11678

For research use only

Overview

Product Name	ADCY3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	125kDa
GeneID?Human?	ADCY3 KIAA0511
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>adenylate cyclase 3(ADCY3) Homo sapiens This gene encodes adenylyl cyclase 3 which is a membrane-associated enzyme and catalyzes the formation of the secondary messenger cyclic adenosine monophosphate (cAMP). This protein appears to be widely expressed in various human tissues and may be involved in a number of physiological and pathophysiological metabolic processes. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016],</p>