

RAD9B Polyclonal Antibody

Cat No: HR1AP11675

For research use only

Overview

Product Name	RAD9B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	RAD9B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the rad9 family.,subunit:Interacts with HUS1, HUS1B, RAD1, RAD9A and RAD17.,tissue specificity:Expressed in testis and skeletal muscle.,