

# SIVA Polyclonal Antibody

Cat No: HR1AP11676

For research use only

## Overview

Product Name	SIVA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	SIVA1 SIVA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SIVA1 apoptosis inducing factor(SIVA1) Homo sapiens This gene encodes a protein with an important role in the apoptotic (programmed cell death) pathway induced by the CD27 antigen, a member of the tumor necrosis factor receptor (TNFR) superfamily. The CD27 antigen cytoplasmic tail binds to the N-terminus of this protein. Two alternatively spliced transcript variants encoding distinct proteins have been described. [provided by RefSeq, Jul 2008],</p>